Proofed by: _____ Date: _____

Dimensions checked: _____ Copy checked: _____

Issue Date:

Date:

DCN No:

Signature:

Release to:

Artwork No: CF36-1832

Business Unit: V. Mueller® / Snowden-Pencer Products

Label: Instructions for use

Description: Snowden Pencer Ring Handle (Diamond-Line I), Electrosurgical 2vV, Non Electrosurgical

1.0 ▲ DESIGNATED INSPECTION DIMENSION OR CHARACTERISTIC

2.0 PHYSICAL REQUIREMENTS:

2.1 Dimensions: 4.0" W x 10.0" H

2.2 Material: #40 Opaque

2.3 Color of Ink: BLACK

3.0 SUPPLIER IDENTIFICATION REQUIREMENTS:

Identify each carton or bundle in a shipment with CareFusion stock number,

CareFusion P.O. number, manufacturer's name, and quantity.

4.0 COPY REQUIREMENTS:

4.1 Printed content must match artwork exactly. ▲ See artwork below.

4.2 Units must be free of dirt and smeared ink ▲ which is easily visible at a distance of 12 - 18".

4.3 Clean sharp edges, free of shreds, wrinkles, and loose particles.

4.4 Printing must be legible. ▲

5.0 SPECIAL REQUIREMENTS: Booklet format

Snowden-Pencer™ Endoscopic Products

en	Diamond-Line™ Endoscopic Products
fr	Diamond-Line Produits endoscopiques
de	Diamond-Line Endoskopieprodukte
it	Diamond-Line Prodotti endoscopici
es	Diamond-Line Productos para endoscopia
pt	Diamond-Line Produtos endoscópicos
nl	Diamond-Line endoscopische producten
sv	Diamond-Line Endoskopiska produkter
da	Diamond-Line endoskopiske produkter
fi	Diamond-Line Endoskooppiset tuotteet
no	Diamond-Line Endoskopiprodukter
el	Diamond-Line Προϊόντα ενδοσκόπησης
ja	Diamond-Line 内視鏡製品
cs	Diamond-Line Endoscopické produkty
pl	Diamond-Line Produkty do endoskopii
tr	Diamond-Line Bedenici Bakı Ürünleri
ro	Diamond-Line Produse endoscopice

USA Rx Only

Snowden-Pencer and Diamond-Line are trademarks or registered trademarks of CareFusion Corporation, or one of its subsidiaries.

© Copyright 2012, CareFusion Corporation, or one of its subsidiaries. All rights reserved.

CareFusion
1500 Waukegan Road
McGaw Park, IL 60085 U.S.A.

CF36-1832 / 26-2900-A • 2012-05

CareFusion

Proofed by: _____ Date: _____


Dimensions checked: _____ Copy checked: _____

en

Catalog Numbers

All products covered by these instructions are listed in the appendix.

These instructions for use are intended to be used with the ring handle (Pictured Below)

Indications For Use

Endoscopic devices are designed to transmit cutting, grasping and dissecting force through delicate mechanisms to delicate working tips in minimally invasive cardiac, general surgical, and plastic surgery.

How Supplied

Snowden-Pencer devices are packaged as non-sterile. Cleaning and sterilization must occur prior to use.

Limitations on Reprocessing

Repeated reprocessing has minimal effect on these devices. End of life is normally determined by wear and damage due to use.

Warnings

Devices shall be used in accordance with these instructions for use. Read all sections of this insert prior to use. Improper use of this device may cause serious injury. In addition, improper care and maintenance of the device may render the device non-sterile prior to patient use and cause a serious injury to the patient or health care provider.

When not in use, store in a location that is isolated from the patient.

Prior to use of the device in patients who have cardiac pacemakers or other active implants, a possible hazard exists because interference with the action of the device or damage to the device may occur. Consult device manufacturer's instructions for use. In case of doubt, qualified advice should be obtained.

Gas embolism may result from over-insufflation of air, inert gas prior to HF surgery or laser assist gas.

The Snowden Pencer instrument is a monopolar device and should only be connected to a generator that is compatible with monopolar devices.

Do not exceed the 0.625 kVp maximum rating specified for this device. Adjust the electrosurgical generator to a setting appropriate for the maximum peak output voltage.

A dispersive electrode, such as a grounding pad, should be used with the device and generator to prevent burns/injury to the patient when using the monopolar device. Ensure the pad is properly sized and properly placed and always recheck the pad when patient is moved. It is recommended that the device is used with an HF generator that contains contact quality monitoring (return electrode monitoring) with a signal to indicate there is contact to the patient.

Conductive fluids, (e.g., blood or saline) in direct contact with an active electrode may carry electrical current or heat, which may cause unintended burns to the patient.

Do not use Snowden Pencer devices with energy in the presence of combustible/explosive gases or other flammable anesthetics, such as nitrous oxide (N₂O) and oxygen. Non-flammable agents used for cleaning and disinfecting, or as solvents of adhesive should be allowed to evaporate before the application of HF surgery. Attention should be called to the danger of ignition of endogenous gases. Some materials, for example cotton, wool, and gauze, when saturated with oxygen may be ignited by sparks produced in normal use of the HF surgical equipment.

There is a risk of pooling of flammable solutions under the patient or in body depressions such as the umbilicus, and in body cavities such as the vagina. Any fluid pooled in these areas should be mopped up before HF surgical equipment is used.

To reduce capacitive coupling, the device should only be activated when in position to deliver energy to the target tissue.

Activating the electrosurgical unit while simultaneously utilizing the suction/irrigation may alter the path of the electrical energy away from target tissue.

Improper connections of accessories may result in inadvertent accessory activation or other potentially hazardous conditions.

Connect adaptors and accessories to the electrosurgical unit only when the unit is off. Failure to do so may result in an injury or electrical shock to the patient or operating room personnel.

Start with the lowest possible power setting on the electrosurgical generator. To achieve the desired cutting and coagulation, check the patient circuit and then gradually increase the power setting. **Note:** the output power selected should be as low as possible for the intended purpose.

Do not activate the electrosurgical unit until the device has made contact with the patient.

Interference produced by the operation of the HF surgical equipment may adversely influence the operation of other electronic equipment.

Use of low frequency power could result in neuromuscular stimulation.

Proofed by: _____ Date: _____

Dimensions checked: _____ Copy checked: _____

Cautions

Always use the device in field of view of the camera and use caution when using other devices so as not to damage device insulation.

Smoke-plume extraction may be necessary during electrosurgical cutting or coagulation. Device (including insulation) and device accessories must be inspected prior to use to ensure integrity. In particular, electrode cables and endoscopically used accessories should be checked for possible damage to the insulation.

To prevent the possibility of electrical shocks or burns, do not use devices with breaks in the insulation.

When used in conjunction with laser equipment, always use caution (i.e. eye protection).

If there are any variations between these instructions and either your facility's policies and/or your cleaning/sterilization equipment manufacturer's instructions, those variations should be brought to the attention of the appropriate responsible hospital personnel for resolution before proceeding with cleaning and sterilizing your devices.

Use of device for a task other than that for which it is intended will usually result in a damaged or broken device.

Examples

1. Use of a delicate dissector as a grasper.
2. Use of a delicate scissor to cut suture.
3. Use of a dissector to remove clips.
4. Use of a 5mm grasper or dissector instead of a 10mm claw extractor forceps to remove excised tissue through cannula.

Prior to use, inspect device to ensure proper function, insulation and condition. Do not use devices if they do not satisfactorily perform their intended function or have physical damage.

Inspect insulation. Any interruptions in the coating may compromise the safety of the device. To prevent the possibility of electrical shocks or burns, do not use devices with breaks in the insulation.

Avoid mechanical shock or overstressing the devices. Close distal ends prior to insertion or removal through cannulas.

Devices are designed to be held with one finger and thumb in ring handles. If force is applied in any other fashion (such as holding handles in a whole hand pistol grip), the mechanism can be easily overstressed resulting in damage or breakage.

Only the cleaning and sterilization processes which are defined within these instructions for use have been validated.

Always use caution when inserting or removing devices through cannula. Lateral pressure on the device during removal can damage the working tip, shaft of the device and/or insulation. Be sure the tips are closed and the device is pulled straight out until completely clear of cannula to avoid catching the valve assemblies in cannulas or dislodging the cannula.

Use only neutral pH (6-8) detergent solutions.

Rated Voltage

The maximum recurring peak voltage for electrosurgical devices applicable to this IFU is 0.625 kVp.

Note: Failure of the electrosurgical generator may result in an unintended increase in output power.

Electrode Cable Compatibility

All of Snowden-Pencer monopolar electrosurgical devices are compatible with the Unipolar Bovie Cord with universal plug (88-9199) or equivalent.

Pre-processing Instructions

Initiate cleaning of device within 2 hours of use.

Transport devices via the institutions established transport procedure.

Remove excess gross soil as soon as possible after use by rinsing or wiping the device.

All devices must be processed in the completely open configuration.

All flush ports shall remain in the fully open position.

Manual Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Prepare the enzymatic / neutral pH detergent solution, utilizing tap water with a temperature range of 27°C to 44°C (81°F to 111°F), per manufacturer's instructions.
3. Place device in the open/relaxed position, with flush port open, and completely immerse in the detergent solution and allow device to soak for a minimum of 5 minutes. Actuate all movable parts during the initiation of the soak time.
4. Using a soft bristled brush, remove all visible soil from the device. Actuate device while brushing, paying particular attention to hinges, crevices and other difficult to clean areas. **Note:** It is recommended that the detergent solution is changed when it becomes grossly contaminated (bloody and/or turbid).
5. For lumen devices, use a soft bristled brush with a brush diameter and length that is equivalent to lumen diameter and length. Scrub the lumen (i.e. angulated/nonangulated positions) until no visible soil is detected in the lumen rinsing step below.

Proofed by: _____ Date: _____

Dimensions checked: _____ Copy checked: _____

MA Gogots 1/11/12

6. For lumen devices, place the device into the open/relaxed position with the distal tip pointed down. Flush the device with a minimum of 50 ml of detergent solution utilizing a temperature range of 27°C to 44°C (81°F to 111°F), by using the flushing port located on the handle/shaft. Repeat the flush process a minimum of 2 times (i.e. total of 3 times), ensuring all fluid exiting the lumen is clear of soil.
7. For lumen devices, if visible soil is detected during the final lumen flush, re-perform brushing and flushing of the lumen.
8. Rinse the device by completely immersing in tap water with a temperature range of 27°C to 44°C (81°F to 111°F), for a minimum of 30 seconds to remove any residual detergent or debris.
9. For lumen devices, following the rinsing step, place the device into the open/relaxed position with the distal tip pointed down. Flush the device with a minimum of 50 ml of tap water utilizing a temperature range of 27°C to 44°C (81°F to 111°F), by using the flushing port located on the handle/shaft. Repeat the flush process a minimum of 2 times (i.e. total of 3 times).
10. Dry the device with a clean, lint-free towel.
11. For lumen devices, manipulate the device to allow rinse water to drain from the lumen.
12. Visually examine each device for cleanliness.
13. If visible soil remains, repeat cleaning procedure.

Automatic Cleaning

1. Ensure all pre-processing instructions are followed prior to cleaning.
2. Clean the devices via the Automatic Cleaning parameters below.

PHASE	MINIMUM RECIRCULATION TIME	WATER TEMPERATURE	DETERGENT TYPE AND CONCENTRATION (IF APPLICABLE)
Pre-Wash 1	15 Seconds	Cold Tap Water 1°C - 16°C (33°F - 60°F)	N/A
Enzyme Wash	1 Minute	Hot Tap Water 43°C - 82°C (110°F - 179°F)	<ul style="list-style-type: none"> Detergent: pH-neutral/enzymatic detergent Concentration: Per the detergent manufacturer's recommendations
Wash 1	2 Minute	Tap Water 43°C - 82°C (110°F - 179°F)	<ul style="list-style-type: none"> Detergent: pH-neutral cleanser Concentration: Per the detergent manufacturer's recommendations
Rinse 1	15 Seconds	Tap Water 43°C - 82°C (110°F - 179°F)	N/A
Pure Rinse	10 Seconds	Purified Water 43°C - 82°C (110°F - 179°F)	N/A
Drying	N/A	N/A	N/A

3. Manipulate the device to allow rinse water to drain from the lumen.
4. If visible moisture is present dry the instrument with a clean, lint-free towel.
5. Visually examine each instrument for cleanliness.
6. If visible soil remains, repeat cleaning procedure.

Inspection/Maintenance

Proper care and handling is essential for satisfactory performance of any surgical device. The previous cautions should be taken to ensure long and trouble-free service from all your surgical devices. Inspect devices before each use for broken, cracked, tarnished surfaces, movement of hinges, and chipped or worn parts. If any of these conditions appear, do not use the device. Return devices to an authorized repair service center for repair or replacement.

Before sterilizing, lubricate the device with instrument milk or a steam permeable/ water soluble lubricant, following the lubricant manufacturer's instructions.

Immediately proceed to recommended steam sterilization. Do not store instrument prior to sterilization.

Let devices drip dry for three (3) minutes before packaging for sterilization.

Packaging

Devices can be loaded into dedicated packaging systems. Sterilization wrap material must be cleared for the applicable sterilization modality by your country's regulatory body. Use in accordance with packaging manufacturer's sterilization instructions being sure to protect jaws and cutting edges from damage.

Sterilization

All devices must be processed in the completely open position (i.e. jaws, flush-ports etc). Devices should not be disassembled.

All flush ports shall remain in the fully open position.

All devices shall be positioned to allow sterilant contact of all surfaces.

All devices with concave surfaces shall be configured so that water pooling does not occur.

Proofed by: _____ Date: _____
 Dimensions checked: _____ Copy checked: _____

MA Gogots 5/24/12

Prevacuum Steam Sterilization Parameters

Minimum Preconditioning Pulses: 3
 Minimum Temperature: 132°C (270°F)
 Minimum Exposure Time: 3 minutes
 Minimum Dry Time: 30 minutes
 Sterilization Configuration: Wrapped (2 layer 1-ply or 1 layer 2-ply)

Gravity Steam Sterilization Parameters

Minimum Temperature: 132°C (270°F)
 Minimum Exposure Time: 15 minutes
 Minimum Dry Time (lumens): 45 minutes
 Minimum Dry Time (non-lumens): 30 minutes
 Sterilization Configuration: Wrapped (2 layer 1-ply or 1 layer 2-ply)

Storage

After sterilization, devices must remain in sterilization packaging and be stored in a clean, dry environment.

Warranty

CareFusion guarantees every surgical device bearing the Snowden-Pencer brand name to be free of functional defects in workmanship and materials when used normally for its intended surgical purpose. Any Snowden-Pencer device proving to be defective will be replaced or repaired at no charge. This device carries a lifetime warranty against manufacturer defects, 3 year warranty against wear.

Repair Service

Regardless of age, if any Snowden-Pencer device needs service, return it to an authorized repair service center. For repairs outside the U.S., please contact your local distributor.

Note: All devices being returned for maintenance, repair, etc. must be cleaned and sterilized per these instructions for use prior to shipment.

IEC Classification: BF Applied Part

Contact Information

CareFusion
 1500 Waukegan Rd
 McGaw Park, IL 60085 U.S.A.
 800-323-9088

www.carefusion.com

For domestic inquiries email: GMB-VMueller-Cust-Support@carefusion.com

For international inquiries email: GMB-SIT-International-Team@carefusion.com

Other Resources: To learn more about sterilization practices and what is required of manufacturers and end users, visit www.aami.org, www.aorn.org or www.iso.org

Appendix

Table 1 Non-Electrosurgical Devices 						
All product codes covered in these instructions are listed in the following table.						
89-2343	89-2351	89-2353	89-2355	89-2357	89-2359	89-2371
89-2373	89-2375	89-2377	89-2379			
SP90-2000	SP90-2001	SP90-2002	SP90-2003	SP90-2004	SP90-2005	SP90-2006
SP90-2007	SP90-2008	SP90-2100	SP90-2101	SP90-2102	SP90-2103	SP90-2104
SP90-2106	SP90-2107	SP90-2108	SP90-2200	SP90-2201	SP90-2202	SP90-2203
SP90-2204	SP90-2206	SP90-2207	SP90-2208	SP90-3000	SP90-3001	SP90-3002
SP90-3003	SP90-3004	SP90-3005	SP90-3006	SP90-3007	SP90-3008	SP90-3009
SP90-3010	SP90-3011	SP90-3012	SP90-3013	SP90-3014	SP90-3015	SP90-3016
SP90-3017	SP90-3018	SP90-3019	SP90-3020	SP90-3021	SP90-3022	SP90-3023
SP90-3024	SP90-3026	SP90-3027	SP90-3099	SP90-3101	SP90-3102	SP90-3103
SP90-3104	SP90-3105	SP90-3106	SP90-3107	SP90-3108	SP90-3109	SP90-3110
SP90-3111	SP90-3112	SP90-3113	SP90-3114	SP90-3115	SP90-3116	SP90-3117
SP90-3118	SP90-3120	SP90-3122	SP90-3124	SP90-3200	SP90-3201	SP90-3202
SP90-3203	SP90-3204	SP90-3205	SP90-3206	SP90-3207	SP90-3208	SP90-3209
SP90-3210	SP90-3211	SP90-3212	SP90-3213	SP90-3214	SP90-3215	SP90-3216
SP90-3217	SP90-3218	SP90-3219	SP90-3220	SP90-3221	SP90-3222	SP90-3223
SP90-3224	SP90-3226	SP90-3227				
			CareFusion France 309 S.A.S. 8 bis rue de la Renaissance 44110 Châteaubriant - France		 CareFusion 1500 Waukegan Road McGaw Park, IL 60085 U.S.A.	

Table 2: Electrosurgical Devices - Devices not CE marked

89-2361	89-2381	SP90-4000	SP90-4001	SP90-4002	SP90-4003	SP90-4004
SP90-4005	SP90-4006	SP90-4007	SP90-4008	SP90-4009	SP90-4010	SP90-4011
SP90-4012	SP90-4013	SP90-4014	SP90-4100	SP90-4101	SP90-4102	SP90-4103
SP90-4104	SP90-4105	SP90-4110	SP90-4112	SP90-4113	SP90-4114	SP90-4200
SP90-4201	SP90-4202	SP90-4203	SP90-4204	SP90-4205	SP90-4206	SP90-4207
SP90-4208	SP90-4209	SP90-4210	SP90-4211	SP90-4212	SP90-4213	SP90-4214
SP90-5004	SP90-5005	SP90-5105	SP90-5204	SP90-5205		

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